



medical **ECONET**  
GERMANY



## ECOSound

Easy-to-use foetal pocket Doppler

- Large LCD display for showing heart rate and battery status
- Highly sensitive 2 MHz probe, waterproof
- Integrated loudspeaker
- Audio output for headphones
- Noise reduction for optimum signal quality
- Approx. 6 hours of operation with continuous measurement
- Automatic switch-off

# Specifications

Model	ECOsound
Display	LCD display with battery indicator
Heart rate measurement range	50-240 bpm
Accuracy	± 2 %
US frequency	2 MHz
US intensity	<10mW/cm2
Loudspeaker	1.2 W integrated loudspeaker
Audio output for headphones	yes
Automatic switch-off	after 5 minutes
Power supply	2 x AA batteries, 1.5V
Sensitivity	from the 10th week of pregnancy

