



ECOsound

Easy-to-use pocket size fetal doppler

- Large LCD for heart rate and battery status
- Ultra sensitive 2 MHz transducer
- Built in loud speaker
- Audio output-jack for earphone
- Active noise reduction to optimize signal quality
- Auto shut off to prolong battery life (AA 1.5Vx2)



medical **ECONET**
GERMANY

Specifications

Model	ECOsound
Display	LCD battery display Signal quality display
Heart rate range	50-240 bpm
Accuracy	± 2 %
Ultrasound frequency	2 MHz
Ultrasound intensity	<10mW/cm2
Speaker output power	1.2 W Built-in loud speaker
Audio output for earphone	0
Auto shut-off	5 minutes
Power	1.5V battery x 2 (type: AA)
Sensitivity	10 weeks onward



medical **ECONET**

Im Erlengrund 20, 46149 Oberhausen - Germany

t: +49 (0)208 / 377 890 - 0 info@medical-econet.com
f: +49 (0)208 / 377 890 - 55 www.medical-econet.com

distributor

*All specifications are subject to change without notice.