



medical **ECONET**  
GERMANY



## ME 7

The efficient lightweight among the ME 7 finger pulse oximeters is ideal for fast and uncomplicated detection of SpO2 and pulse.

### The perfect finger pulse oximeter...

... for measuring oxygen saturation in arterial blood (partial oxygen saturation or SpO<sub>2</sub>).

It shows:

SpO<sub>2</sub> plethysmogram, 2 SpO<sub>2</sub> value, 3 pulse rate value, battery power and pulse intensity.

The oximeter uses two 1.5 V AAA batteries.

The unit does not have an alarm function.

