



medical **ECONET**
GERMANY



mex+ DG70

Dental X-ray generator

- **Reliable:** Latest technology for a reliable application
- **Efficient:** Compact and easy to use, for a fast and certain diagnosis
- **Usability:** Intuitive usable control panel with LED display
- **Performance:** High image quality based on the small focal spot (0,4 mm) and 70 kV
- **Compact:** Space-saving wall-mount version for small rooms
- **Mobility:** Simple movement due to integrated mobile stand


Performance meets mobility.

With its power (70kV, 3.5 mA) and its small focal spot (0.4 mm), the meX+ DG70 dental x-ray unit creates images with excellent resolution, high contrast and brilliant image sharpness.

The user-friendly control panel with LED display simplifies the daily workflow.

Specifications

Model	meX+DG70
Electronic properties	Power supply: 230 V \pm 10% Frequency: 50 Hz
Tube head	Total filtration: \geq 2,5 mm Al. eq Radiation field: @200 \varnothing 55 mm
X-ray tube	Focal spot: 0,4 mm Anode material: tungsten
Radiological properties	High voltage: 70 kVp Anode current: 3,5 mA Exposure time: 0,04 - 2,00



Keypad foils



Vet-Version



Human-Version

medical **ECONET**

Im Erlengrund 20, 46149 Oberhausen - Germany

t: +49 (0) 208 377 890 - 0 info@medical-econet.com
f: +49 (0) 208 377 890 - 55 www.medical-econet.com
facebook.com/medicaleconetgermany
twitter.com/medicaleconet

CE 0051

Dealer

*All specifications are subject to change without notice.